Pectoral Implants (including PowerFlex™ Implants)

PHYSICIAN'S INFORMATION DATA SHEET

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SPA

POR

Canadian Physicians, Please go to www.implantech.com/Canada to view the IFUs for Canadian Physicians.

Symbol Definition - Refer to label for symbols applicable to your device.

	Caution, Consult Instructions for Use
	Do Not Reuse, "Single Use Only"
	Catalog Number
	Batch Code, "Lot"
	Use By

	Sterilized Using Steam or Dry Heat
	Non-Sterile
	Manufactured By
	Manufacture Date
	U.S. Federal law restricts this device to sale by or on the order of a physician.

PHYSICIAN'S INFORMATION DATA SHEET

DESCRIPTION

Implantech's Pectoral Implants, including PowerFlex™ Implants and Sizers are 100% medical grade solid silicone elastomers that have been molded into various convex and oval shapes and durometers. These long-term implants are designed as space-occupying devices for augmentation and reconstructive surgery.

INDICATIONS FOR USE

Pectoral Implants, including PowerFlex™ are indicated for augmentation or reconstruction of the pectoralis region.

Sizers function as space-occupying devices for temporary insertion intraoperatively to approximate the shape and size of the associated implants. They are colored blue to show that they are not intended for long-term implantation.

CONTRAINDICATIONS

1. Lack of sufficient tissue covering
2. Autoimmune disease such as Lupus Erythematosus or Rheumatoid Arthritis
3. Infectious disease
4. Diabetes or other endocrine disease, which compromises circulation
5. The presence of any medical or psychiatric condition that would cause a poor risk for this type of surgery. The surgeon must determine whether the patient is suitable for this type of surgery.

NOTE: EVERY PATIENT SHOULD BE MADE AWARE OF THE FOLLOWING FACTS.

Augmentation surgery can provide satisfaction to patients. Nevertheless, it is not without potential complications and risks. Implantation surgery is an **ELECTIVE PROCEDURE** and the patients should be well counseled on the risk-benefit relationship. The possibility of explant surgery taking place at any time after implantation should be discussed with the patient.

IMPLANTECH RELIES ON THE SURGEON TO ADVISE THE PATIENT OF THE COMPLICATIONS AND RISK ASSOCIATED WITH BOTH THE IMPLANT AND THE SURGERY ITSELF.

COMPLICATIONS AND WARNINGS

Resorption of bone underlying the implant can occur.

Inadequate tissue covering, inadequate size surgical pocket or too large an implant can result in tissue necrosis and extrusion of the implant.

Displacement or shifting of the implant can occur from too large a pocket.

Formation of a fibrous tissue capsule surrounding the implant is a normal physiologic response to any foreign body. Contracture of the fibrous tissue capsule can result in firmness, implant displacement or buckling, discomfort, or pain.

Utmost care should be taken to avoid damage to the implant during handling or in surgery. Care should be used when trimming to avoid creating cuts or tears into the device which may propagate (trimmings should be continuous and clean cut). Avoid contact with sharp objects such as surgical instruments or suture needles that may produce cuts or tears into the device. Avoid undue handling with blunt instruments or manipulation which can tear the device.

Complications associated with all surgical procedures should be discussed with the patient (e.g. infection; poor reaction to medication/surgical procedures; poor wound healing; hematoma; serous fluid accumulation; nerve damage or irritation; neuralgia; loss of sensation; patient intolerance to any foreign implant).

Surgeons should consider the following prior to any surgical implant procedure:

1. Known allergies to silicone materials.
2. Active infection.
3. Compromised circulation to the skin and soft tissues in areas to be implanted.
4. Severe medical problems relating to the areas to be implanted.
5. Unrealistic patient expectations.

Due to the wide variety of patients' physical responses to implant surgery and the differences in surgical techniques and medical treatments, as well as the possibility of complications or trauma, **PATIENTS SHOULD BE ADVISED THAT THESE SHOULD NOT BE CONSIDERED LIFETIME IMPLANTS AND EXPLANT SURGERY MAY BE INDICATED AT ANY TIME. IMPLANTECH MAKES NO REPRESENTATIONS FOR THE TERM OF IMPLANTATION OF THE DEVICE.**

Implants and Sizers are intended for single patient use and single use only.

WARNING: UNDER NO CONDITIONS SHOULD AN EXPLANTED PRODUCT BE REIMPLANTED!

This is because recommended recleaning and resterilization may not adequately remove biological residues such as blood, tissue, and other matter, which could retain resistant pathogens.

SURGICAL PROCEDURE

Proper surgical procedures and techniques are the responsibility of the medical profession. Each surgeon must evaluate the suitability of the procedure based upon current accepted techniques, individual judgement, and experience.

Proper size and shape of implants must be determined for the individual patient by the surgeon. Sizers are available from Implants for the convenience of the surgeon.

SUPPLIED STERILE OR NON-STERILE

Implants supplied in sterile form are processed by a validated, strictly controlled cycle. Sterility is verified in accordance with AAMI / ANSI / ISO Standards. Sterility of the implant is maintained only if the package is intact and undamaged.

This device is sterilized by a Dry Heat process which results in a natural variation in translucency from product to product.

If the sterile implant becomes contaminated prior to implantation, it may be recleaned and resterilized by autoclaving according to procedures described below. The implant should not be resterilized by ethylene oxide gas, since excessive accumulation of ethylene oxide residuals may cause adverse tissue reaction.

Implants supplied in non-sterile form must be sterilized by autoclaving before use. (SEE STERILIZATION INSTRUCTIONS).

PACKAGING

Sterile product is supplied in a sealed double package. Sterility is not guaranteed if the package has been damaged or opened.

Non-sterile product is supplied in a sealed, single package.

A tear-off patient chart record is provided on both sterile and non-sterile packages.

TO OPEN STERILE PRODUCT:

1. Peel open outer pouch under clean aseptic conditions.
2. Invert outer pouch over sterile fields, allowing sealed inner package to fall on field.
3. Using aseptic precautions, peel open inner package and retrieve implant.

NOTE: ATTACH PATIENT RECORD PORTION OF THE LABEL TO PATIENT'S CHART.

NON-STERILE PRODUCT AND SIZERS

1. If sterilization or resterilization is required, the following cleaning and sterilization techniques have been found effective.
2. Remove implant from its package in a clean environment using gloved hands. Do not sterilize in the packaging system supplied. Utmost caution should be taken to avoid contamination as lint, fingerprints, talc, and other surface contaminants can cause foreign body reactions.
3. **TO CLEAN:** Wash thoroughly in a hot water-soap solution. Use a mild non-oily soap. Rinse thoroughly in hot water followed by distilled water and sterilize.

DO NO USE SYNTHETIC DETERGENTS!!

DO NOT USE OIL BASED SOAPS!!

4. Sterilize per **STERILIZATION INSTRUCTIONS** below.

NOTE: ATTACH PATIENT RECORD PORTION OF LABEL TO PATIENT'S CHART.

STERILIZATION INSTRUCTIONS

The Sterilization Instructions provided are to be used only as a guide. We recommend that the efficacy of the specific autoclaving cycle employed be established by appropriate methods such as the use of commercially available chemical or biological monitors.

DO NOT USE ETHYLENE OXIDE STERILIZATION.

Wrap implant in a clean, lint-free material and place in a clean open autoclaving tray. Sterilize by the following gravity displacement autoclaving cycle.

1. Minimum of 30 minutes at 270°F, 30 psi.

DO NOT USE AN AUTOCLAVE CYCLE THAT INCLUDES A VACUUM PHASE (SUCH AS PRE-CONDITIONING PHASE OR IN POST-CONDITIONING AS A DRYING PHASE).

When steam autoclaving, the user should verify that the sterility procedure will produce a sterility assurance level of 10^{-6} by using an appropriate biological indicator and validation program.

WARNING: IMPLANTS SHOULD BE ALLOWED TO COOL THOROUGHLY AFTER AUTOCLAVING.

PRODUCT EXAMINATION AND HANDLING

1. Prior to implantation, products should be visually examined for any evidence of particulate contamination or damage. **DAMAGED PRODUCTS SHOULD NOT BE IMPLANTED.** Do not attempt to repair damaged products.
2. The implant should be kept submerged in sterile normal saline prior to implantation to prevent contact with airborne and surgical field particulate contaminants.
3. Care must be taken to prevent possible surface contamination by talc, dust and skin oils which might adversely affect the suitability of the implant.

WARRANTY, LIMITATION OF LIABILITY AND DISCLAIMER OF OTHER WARRANTIES.

Implantech warrants that reasonable care was used in the manufacture and production of the product. Because Implantech has no control over the conditions of use, patient selection, surgical procedure, post-surgical stresses, or handling of the device after it leaves our possession, Implantech does not warrant either a good effect or against an ill effect following its use. Implantech shall not be responsible for any incidental or consequential loss, damage, or expenses directly or indirectly arising from use of this product. Implantech's **SOLE** responsibility, in the event Implantech determines the product was defective when shipped by Implantech shall be replacement of the product. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for use.

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© 2024 Implantech
Implantech Associates, Inc.