Silicone Sheeting

Physician's Information Data Sheet

SILICONE SHEETING

PHYSICIAN'S INFORMATION DATA SHEET

English

PRODUCT DESCRIPTION:

Implantech's AlliedSil™ Silicone Sheeting are sheets of vulcanized silicone ranging from 0.127mm (0.005") to 1.524mm (0.060") in thickness, and are available in sheet sizes ranging from 50.8mm x 76.2mm (2" x 3") to 304.8mm x 431.8mm (12" x 17"). Most sheets are available either non-reinforced or reinforced with a thin polyester mesh material. Sheeting labeled as short term is restricted to less than a 30 day term of implantation, and sheeting labeled as long term is suitable for an indefinite term of implantation. (NOTE: All configurations of AlliedSil™ Silicone Sheeting are suitable for use in the Magnetic Resonance (MR) environment.)

Silicone Sheeting, Short Term (Reinforced and Non-Reinforced)

(Note: Short Term Sheeting is limited to less than 30 days.)

INTENDED USE:

The intended use for this silicone sheeting is as a short term space-occupying substance for augmentation or reconstructive surgeries of the head and neck, for management of keloid or hypertrophic scars, or as a wound or tissue covering.

Silicone Sheeting, Long Term (Reinforced and Non-Reinforced)

INTENDED USE:

The intended use for this silicone sheeting is as a space-occupying substance for augmentation or reconstructive surgeries of the head and neck.

CONTRAINDICATIONS

Lack of sufficient tissue covering.

This product was designed for use by licensed physicians with proper training and experience. Implantech relies on the user to carefully review the medical literature of indications, contraindications and risks in using this product. It is the responsibility of the using physician to consider these risks before prescribing this product, and to advise the patient of the complications and risk associated with the implant and the surgery.

Implants are intended for **single patient use only**. **WARNING: UNDER NO CONDITIONS SHOULD AN EXPLANTED PRODUCT BE REIMPLANTED!** Re-use following recleaning and resterilization may not adequately remove biological residues such as blood, tissue, and other matter, which may retain resistant pathogens.

COMPLICATIONS AND WARNINGS:

1. The surgeon should be aware of the complications associated with implantation of silicone devices. Risks include formation of fibrous tissue capsule around the implant, cosmetic deformation, displacement or shifting, bone resorption underlying the implant, extrusion and common surgical risks including pain, infection, hematoma, serous fluid accumulation, nerve damage or irritation, neuralgia, loss of sensation, poor wound healing, malaise, fever, foreign body reaction or patient intolerance to a foreign implant, and repeated surgery.
2. Utmost care should be taken to avoid damage to the implant during handling or in surgery, particularly with the use of sharp objects. Avoid undue handling with blunt instruments or manipulation.
3. It is the responsibility of the physician to advise patients of risks, complications and benefits of implant surgery.
4. Failure to properly inspect implant and determine it is suitable prior to surgery, may result in untimely cancellation of surgery.

(NOTE: There are no special requirements for disposal of unused devices. Techniques for removal of devices are dependent on the particular intended use and are the responsibility of the surgeon. Disposal of explanted devices should be handled in accordance with accepted protocols for handling potentially bio-hazardous materials.)

INSTRUCTIONS FOR USE:

Proper surgical procedures and techniques, including the determination of whether and how to fixate the device, are the responsibility of the medical profession.

INSPECTION UPON RECEIPT:

1. Make sure package is unopened and undamaged. Do not use if damaged.
2. Read labels to ensure the package contains the required product.
3. Carefully check size and style through package window.

PRESTERILIZED PRODUCTS:

1. Products labeled "Sterile" have been gamma radiation sterilized in a double Tyvek pouch.
2. Confirm package integrity and sterility have not been compromised.
3. Just prior to surgery, carefully remove outer Tyvek pouch and aseptically transfer sterile inner pouch to sterile field.
4. Using accepted sterile surgery techniques, peel pouch and place implant into a basin of sterile saline to protect from contamination.

5. Silicone implants are made of medical implant grade quality silicone known for excellent biocompatibility. However, contaminants such as dust, talc, detergents, antiseptic agents and skin oils transferred to the product during handling can evoke a foreign body response.

NON-STERILE PRODUCT:

1. If sterilization or resterilization is required, the following cleaning and sterilization techniques have been found effective.
2. Remove implant from its package in a clean environment using gloved hands. Do not sterilize in the packaging system supplied. Utmost caution should be taken to avoid contamination as lint, fingerprints, talc, and other surface contaminants can cause foreign body reactions.
3. TO CLEAN: Wash thoroughly in a hot water-soap solution. Use a mild non-oily soap. Rinse thoroughly in hot water followed by distilled water and sterilize. **DO NOT USE SYNTHETIC DETERGENTS!! DO NOT USE OIL BASED SOAPS!! DO NOT USE WASHER-DISINFECTOR EQUIPMENT!!**
4. Sterilize per STERILIZATION INSTRUCTIONS below.

NOTE: Attach patient record portion of label to patient's chart.

STERILIZATION INSTRUCTIONS

DO NOT USE ETHYLENE OXIDE STERILIZATION.

Package in an autoclavable pouch or wrap in a clean, lint-free material and place in a clean, open autoclaving tray. Sterilize by one of the following gravity displacement autoclaving cycles:

Minimum of 30 minutes at 132° C (270° F) with a gravity air displacement cycle, OR
Minimum of 4 minutes at 132° C (270° F) with a prevacuum conditioning cycle.

When using autoclaves for sterilization, it has to be strictly ensured that the steam used is absolutely free of foreign substances. (Note: Do not exceed cumulative autoclaving cycles of 600 minutes at 134° C.)

Storage: After autoclaving, the autoclaved device should be stored in such a way as to maintain the sterility of the device until ready for use.

WARNING: IMPLANTS SHOULD BE ALLOWED TO COOL THOROUGHLY AFTER AUTOCLAVING.


The instructions provided above have been validated by Implantech as being CAPABLE of preparing the silicone sheet for use. It remains the responsibility of the user (processor) to ensure that the processing as actually performed using the equipment, materials and personnel in the facility achieve the desired results. This normally requires validation and routine monitoring of the process. We recommend user strictly observe the instructions for use given by the manufacturer of the autoclave.

WARRANTY

Implantech warrants that reasonable care was used in the design and manufacture of the product described herein. Implantech's sole responsibility and liability shall be that it will replace product if shown to be defective. Implantech shall not be liable for any other incidental damages. This warranty is in lieu of and excludes express or implied warranties by operation of law or otherwise, including warranties of merchantability or fitness for use.

Symbol Definition

– Refer to label for symbols applicable to your device.

	Caution, Consult Instructions for Use		Sterilized Using Steam or Dry Heat
	Do Not Reuse, "Single Use Only"		Sterilized Using Irradiation
	Catalog Number		Non-Sterile
	Batch Code, "Lot"		Manufactured By
	Use By		Manufacture Date
	U.S. Federal law restricts this device to sale by or on the order of a physician.		Authorized Representative in the European Community

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