



Medical Grade FEP Fluoropolymer Film

Description:

Implantech FEP Fluoropolymer Films are melt adhesives for heat bonding or laminating materials that have high temperature stability. The FEP films can be used in conjunction with all of Implantech's fluoropolymer products. The films are available in 0.5 mil, 1 mil and 2 mil gauges. Implantech FEP films have been qualified for permanent implant applications in accordance with ISO 10993:2003.

Ancillary Services:

- Custom Configurations
- Product Development
- Contract Assembly

Contact Us For More Information:

Implantech
6025 Nicolle St.
Ventura, CA 93003
800-733-0833, ext. 212
fp@implantech.com

Physical Properties:

Property	Typical Result	Method
Appearance	Clear	Visual
Melting Point	265 C	ASTM D-2116
Specific Gravity	2.14 g/cc	ASTM D-2116
Tensile Strength	3600 psi.	ASTM D-2116

