



## Medical Grade Non Reinforced ePTFE Sheets

### Description:

Implantech Non Reinforced ePTFE Sheets are soft and flexible sheets of PTFE with a controlled porous microstructure. The ePTFE sheets are manufactured in an USFDA-approved medical device manufacturing facility with the highest standards of purity. Tightly controlled process specifications and lot testing ensure that each and every product conforms to specification. Implantech ePTFE materials have been qualified for permanent implant applications in accordance with ISO 10993:2003.

### Ancillary Services:

- Custom Configurations
- Product Development

### Contact Us For More Information:

Implantech  
 6025 Nicolle St.  
 Ventura, CA 93003  
 800-733-0833, ext. 212  
[fp@implantech.com](mailto:fp@implantech.com)

### Physical Properties:

Property	Typical Result
Nodal Distance	30 microns
Density	0.35 g/cc
Tensile Strength	1700 psi

