



NEWS FOR YOUR PRACTICE

APRIL 2021

800.733.0833 | implantech.com



When is a filler not an injectable?

PICTURES SPEAK VOLUMES

Check out some before-after photos

TRAINING WORKSHOPS

Sign up for updated class schedules

FANCY MEETING YOU

Meetings where we hope to see you!



Free practice enhancement is just a click away!

Along with our industry-leading aesthetic products, Implantech offers a number of ways to assist you and your practice.

Just create an account as either a **USA** or **International** physician. Once we verify it, you can log in whenever you wish and click on **Resources**.

In addition to other info, you'll have access to some terrific start-to-finish procedure

When is a filler not an injectable?

The short answer? When the filler is a **facial implant** from Implantech.

Along with sagging skin, a loss of soft tissue "fullness" is a typical hallmark of the aging face. Thus, restoring volume is critical for recapturing a more aesthetic, youthful look.

Often, the same dermal fillers used to lessen wrinkles are also injected to create facial volume. However, follow-up injections are needed to replace the filler material eventually resorbed (sometimes unevenly) by the body.

This resorption—along with potential complications such as asymmetry, lumps or palpability under the skin—can make fillers less than ideal for adding rejuvenating volume to large facial areas such as chins, cheeks and jaws.

For a favorable alternative, turn to Implantech.



Conform Terino Malar Shell® Implant
Photos courtesy of Joe Niamtu, III, DMD

Pre-formed and available in a variety of

sizes, **our facial implants** can help to re-


augment and re-position facial tissue lacking its original volume.

For example:



Our **Instructional Videos** cover implant procedures for all areas of the face, including midfacial,

chin, jaw, premaxillary, paranasal and temporal. Also available are videos for body contouring with gluteal, pectoral and calf implants.



Check out the **Surgical Tips & Techniques** section to access helpful info from the doctors who designed some of our facial implants.

Discussed in the technique guides are preparation, incision and placement along with indications and contraindications.

Whether you're a novice or a seasoned pro when it comes to facial and body-contouring implants, you're sure to find invaluable insight within these videos and written materials.

Also available are numerous product information sheets, a complete interactive catalog and useful links.

We add new items on a regular basis. But if you don't find what you're looking for, just contact another important Implantech resource—our friendly customer service team!

The result? A more youthful and aesthetic look for your patient.

Choose from volumizing implants designed to enhance midfacial, chin, jaw and temporal regions. Besides our **flagship silicone line**, we offer **Composite™ implants** made of ePTFE-coated silicone and other implants made of **solid ePTFE**.

Check it all out at implantech.com.

QUESTIONS ABOUT OUR PRODUCTS OR THEIR USES?

implantech.com 800.733.0833

IMPLANTECH | implantech.com



IMP-642/0421