

**Implantech**<sup>®</sup>  
Superior Patient Aesthetics

## NEWS FOR YOUR PRACTICE

NOVEMBER 2022  
800.733.0833 | [implantech.com](http://implantech.com)

# SPOTLIGHT ON

## Mandibular Pre Jowl Chin™ Implant

Augmenting the  
chin and the pre-jowl sulcus

### PICTURES SPEAK VOLUMES

Check out some before-after photos

### TRAINING WORKSHOPS

Sign up for updated class schedules

### FANCY MEETING YOU

Meetings where we hope to see you!

### 3D Accuscan PATIENT-SPECIFIC® IMPLANTS

Cleared by the FDA  
for both face and body

#### **When stock implants just won't do!**

For implants tailor-made for most soft-tissue or bony deficiencies, look no further than our **3D Accuscan® system**.

Simply put, 3D Accuscan Patient Specific® implants can make aesthetic changes to the face and body that just aren't possible with standard implants.

For example, check out the following case and photos, courtesy of Barry L. Eppley, MD, DMD:

#### **The Mandibular Pre Jowl Chin™**

*For augmenting the chin-to-jowl groove AND the chin*

**NOTE:** Last month, we featured our Mandibular Pre Jowl® implant for patients who do NOT require added chin projection. This month's implant is for those that do!

As a result of progressive atrophy of soft tissue in the area between the chin and jowl—and the gradual bone reduction of the mandible in this same area—a patient may develop a groove or *pre-jowl sulcus*.

Though a facelift can tighten this area, the sulcus often remains post-operatively.

Moreover, patients with such a sulcus may also exhibit a weak chin profile.

**The Solution** – Implantech's *Mandibular Pre Jowl Chin™ implant*.




Along with enhanced anterior projection, the MPJC provides a subtle increase in lateral fullness precisely positioned to augment the aging jowl line.



## BEFORE

The patient desired a look with a stronger jawline and a more defined cheekbone appearance.



A 3-D digital model of the affected area was made from the patient's CT scans.

Viewing the model from all angles, the surgeon worked with our technician online in real time to design tailor-made implants.

After approving the design (first in a PDF and later, a sample), the surgeon received the actual solid silicone implant just as envisioned.



## AFTER

The patient's desired look was achieved with patient-specific jawline and infra-orbital-malar implants.

As an option, we can also produce a three-dimensional model of your patient's anatomy as a consultation or surgical aid.

**Contact us today** for more information about the 3D Accuscan Patient-Specific Implant® system.

The resulting augmentation enables a fuller, more youthful and aesthetic look.

### *Glasgold Wafer™ (GLM-2)*

Adds 2mm projection to MPJC



For patients who need it, you can extend the MPJC by 2mm with our Glasgold Wafer™ (cat. no. GLM-2).

As with our Mandibular Pre Jowl® for those not needing additional projection, numerous surgeons note their use of the MPJC with facelifts to maintain the defined lateral jawline following the lift.

### The Placement



As noted is his [2008 article](#), James R. Shire, MD, FACS, states, "Placement of the [implant] is through a submental incision, and is usually performed during the face-lift surgery."

**Take a look** at Dr. Shire's description of his surgical technique, and how critical he believes this implant to be for the aging patient's satisfaction.

### The Result

# STAY TUNED!

We'll soon have 2023 dates for our popular hands-on cadaver workshop!



Advanced Techniques  
Using Facial & Body-Contouring Implants

Questions about our  
products or their uses?

[implantech.com](https://www.implantech.com)

800.733.0833

Before

After



Photos courtesy of Rady Rahban, MD, FACS

Contact [Implantech](https://www.implantech.com) for more information about the Mandibular Pre Jowl Chin™ implant or any of our superior products.

IMPLANTECH | [implantech.com](https://www.implantech.com)



IMP-642/1122