

WELCOME SPECTRUM FACIAL IMPLANT CUSTOMERS

Dear Valued Customer,

As Abbvie has made the decision to discontinue their Spectrum Facial Implant line, they are referring valued customers like you to us—Implantech Associates Inc.—the trusted leader in facial and body contouring implants. We are honored to have the opportunity to support your needs moving forward.

To ensure a smooth and seamless transition, we've created a detailed cross-reference guide to help you identify the recommended Implantech replacements for the products you currently use.

At Implantech, we take pride in manufacturing and delivering premium-quality implants designed for exceptional fit and performance. Whether selecting from our extensive range of stock styles and sizes or opting for a fully customized patient-specific design, Implantech has the right solution for you.


If you have questions or would like personalized assistance, we invite you to complete our quick inquiry form. Once submitted, a dedicated Implantech team member will follow up promptly to provide expert guidance tailored to your needs.

We also encourage you to contact our friendly customer service team directly at CustomerService@implantech.com or by calling 800.733.0833 (option 1).

Thank you for your trust. Whether you're new to Implantech or already a valued customer, we're committed to making this transition as smooth and successful as possible. We look forward to supporting your continued success and are honored to have you as part of the Implantech family.

Warm regards,
Implantech


Anatomical Malar Styles I and II



Spectrum's Anatomical Malar Styles I and II are some of the original designs for malar implants which have served as cornerstones in addressing malar augmentation. Over time, as surgical techniques and anatomical understanding have advanced, surgeons and engineers have worked together to deliver refined and enhanced newer designs.


Catalog Number	Nominal Dimensions		
	A	B	C
S120-214S	4.0cm	2.0cm	0.3cm
S120-224S	4.6cm	2.4cm	0.4cm
S120-234S	5.2cm	2.6cm	0.5cm

- Sizer set S900-010
- Implant Supplied Sterile



Catalog Number	Nominal Dimensions		
	A	B	C
S130-314S	4.7cm	2.5cm	0.5cm
S130-324S	5.0cm	3.2cm	0.6cm
S130-334S	5.6cm	3.4cm	0.7cm

- Sizer set S900-011
- Implant Supplied Sterile



ALTERNATIVE

ALTERNATIVE

ALTERNATIVE



Terino Malar Shell®


The Terino Malar Shell® represents this evolution, crafted with a focus on improving the fit and contour to align more naturally with the zygoma. This updated design reflects the input and expertise of leading surgeons, ensuring that it not only builds on the strengths of its predecessor but also incorporates modern advancements for improved outcomes and precision.

Moderate projection to the malar prominence (3-4mm projection): Implantech Sizes TMS-S3, S4, M3 and M4, L3, L4, XL4


Moderate to high projection to the malar prominence (5-7mm projection): Implantech Sizes TMS-M5, L5, L6, XL5, XL6, XL7

Cat. No.	Size	A	B	C
TMS-S3	Small	4.8	2.3	0.3
TMS-S4	Small	4.8	2.3	0.4
TMS-M3	Medium	5.1	2.8	0.3
TMS-M4	Medium	5.1	2.8	0.4
TMS-M5	Medium	5.1	2.8	0.5
TMS-L3	Large	5.2	3.3	0.3
TMS-L4	Large	5.2	3.3	0.4
TMS-L5	Large	5.2	3.3	0.5
TMS-L6	Large	5.2	3.3	0.6
TMS-XL4	X-Large	6.0	3.7	0.4
TMS-XL5	X-Large	6.0	3.7	0.5
TMS-XL6	X-Large	6.0	3.7	0.6
TMS-XL7	X-Large	6.0	3.7	0.7
TMS-SIZ-NS	Sizes	Set of 0.4cm Projection Sizes		

Dimensions in centimeters




Anatomical Malar Styles I and II



Spectrum's Anatomical Malar Styles I and II are some of the original designs for malar implants which have served as cornerstones in addressing malar augmentation. Over time, as surgical techniques and anatomical understanding have advanced, surgeons and engineers have worked together to deliver refined and enhanced newer designs.

Catalog Number	Nominal Dimensions		
	A	B	C
S120-214S	4.0cm	2.0cm	0.3cm
S120-224S	4.6cm	2.4cm	0.4cm
S120-234S	5.2cm	2.6cm	0.5cm


- Sizer set S900-010
- Implant Supplied Sterile



RIGHT SIDE IMPLANT

Catalog Number	Nominal Dimensions		
	A	B	C
S130-314S	4.7cm	2.5cm	0.5cm
S130-324S	5.0cm	3.2cm	0.6cm
S130-334S	5.6cm	3.4cm	0.7cm

- Sizer set S900-011
- Implant Supplied Sterile



RIGHT SIDE IMPLANT

ALTERNATIVE

ALTERNATIVE


ALTERNATIVE




Conform Terino Malar
Shell®

Cat. No.	Size	A	B	C
CTMS-S4	Small	4.8	2.3	0.4
CTMS-M4	Medium	5.1	2.8	0.4
CTMS-L4	Large	5.2	3.3	0.4
CTMS-XL4	X-Large	5.7	3.4	0.4
CTMS-SIZ-NS	Sizers	Complete Set		

Dimensions in centimeters




Anatomical Malar Styles I and II



Spectrum’s Anatomical Malar Styles I and II are some of the original designs for malar implants which have served as cornerstones in addressing malar augmentation. Over time, as surgical techniques and anatomical understanding have advanced, surgeons and engineers have worked together to deliver refined and enhanced newer designs.


Catalog Number	Nominal Dimensions		
	A	B	C
S120-214S	4.0cm	2.0cm	0.3cm
S120-224S	4.6cm	2.4cm	0.4cm
S120-234S	5.2cm	2.6cm	0.5cm

- Sizer set S900-010
- Implant Supplied Sterile



Catalog Number	Nominal Dimensions		
	A	B	C
S130-314S	4.7cm	2.5cm	0.5cm
S130-324S	5.0cm	3.2cm	0.6cm
S130-334S	5.6cm	3.4cm	0.7cm

- Sizer set S900-011
- Implant Supplied Sterile



ALTERNATIVE

ALTERNATIVE


ALTERNATIVE



ePTFE-Coated Terino
Malar Shell®

Cat. No.	Size	A	B	C
ePTFE-TMS-S	Small	4.8	2.3	0.4
ePTFE-TMS-M	Medium	5.1	2.8	0.4
ePTFE-TMS-L	Large	5.2	3.3	0.4
ePTFE-TMS-XL	X-Large	5.7	3.6	0.4
TMS-SIZ-NS	Sizers	Complete Set (silicone)		

Dimensions in centimeters




Extended Malar



Spectrum’s Extended Malar Implant is a basic oval design intended to add moderate volume to the malar fat pad by only 2mm – 4mm.

Catalog Number	Nominal Dimensions		
	A	B	C
S150-514S	3.6cm	1.8cm	0.2cm
S150-524S	4.0cm	2.0cm	0.2cm
S150-534S	4.3cm	2.3cm	0.4cm

- Sizer set S900-012
- Implant Supplied Sterile



ALTERNATIVE



Terino Malar Shell®


Expertly designed for augmenting patients with malar hypoplasia. With a broad surface area and fine, non-palpable edges, there is a smooth transition between the implant and underlying bone.

Similar projection to the Spectrum Extended Malar (3-4mm projection): Sizes TMS-S3, S4, M3 and M4, L3, L4, XL4.

For a more oval shape, implant can be trimmed to the desired shape.

Cat. No.	Size	A	B	C
TMS-S3	Small	4.8	2.3	0.3
TMS-S4	Small	4.8	2.3	0.4
TMS-M3	Medium	5.1	2.8	0.3
TMS-M4	Medium	5.1	2.8	0.4
TMS-M5	Medium	5.1	2.8	0.5
TMS-L3	Large	5.2	3.3	0.3
TMS-L4	Large	5.2	3.3	0.4
TMS-L5	Large	5.2	3.3	0.5
TMS-L6	Large	5.2	3.3	0.6
TMS-XL4	X-Large	6.0	3.7	0.4
TMS-XL5	X-Large	6.0	3.7	0.5
TMS-XL6	X-Large	6.0	3.7	0.6
TMS-XL7	X-Large	6.0	3.7	0.7
TMS-SIZ-NS	Sizers	Set of 0.4cm Projection Sizers		

Dimensions in centimeters




Profile Malar



Catalog Number	Nominal Dimensions		
	A	B	C
S140-414S	4.6cm	2.3cm	0.4cm
S140-424S	5.1cm	2.7cm	0.5cm
S140-434S	5.4cm	3.2cm	0.5cm
S140-444S	5.7cm	3.5cm	0.5cm

- Sizer set S900-014
- Implant Supplied Sterile



SIMILAR PRODUCT




Terino Malar Shell®

Like Spectrum’s Profile Malar, moderate projection to the malar prominence (4-5mm projection): Implantech Sizes and TMS – S4, M4, L4, XL4, M5, L5, XL5.

Additional sizes available for more or less projection.

Cat. No.	Size	A	B	C
TMS-S3	Small	4.8	2.3	0.3
TMS-S4	Small	4.8	2.3	0.4
TMS-M3	Medium	5.1	2.8	0.3
TMS-M4	Medium	5.1	2.8	0.4
TMS-M5	Medium	5.1	2.8	0.5
TMS-L3	Large	5.2	3.3	0.3
TMS-L4	Large	5.2	3.3	0.4
TMS-L5	Large	5.2	3.3	0.5
TMS-L6	Large	5.2	3.3	0.6
TMS-XL4	X-Large	6.0	3.7	0.4
TMS-XL5	X-Large	6.0	3.7	0.5
TMS-XL6	X-Large	6.0	3.7	0.6
TMS-XL7	X-Large	6.0	3.7	0.7
TMS-SIZ-NS	Sizers	Set of 0.4cm Projection Sizers		

Dimensions in centimeters




Combined Malar



The Combined Malar Implants were designed to provide volume to the midfacial soft tissue with moderate anterior augmentation.

Catalog Number	Nominal Dimensions			
	A	B	C	D
S160-614S	4.5cm	3.0cm	0.4cm	0.35cm
S160-624S	5.0cm	3.3cm	0.45cm	0.4cm
S160-634S	5.6cm	3.7cm	0.5cm	0.45cm

- Sizer set S900-013
- Implant Supplied Sterile



SIMILAR PRODUCT

ALTERNATIVE




Combined Submalar Shell™

Provides more size options than the Spectrum Combined Malar. Augments entire midfacial area by incorporating features of the Terino Malar Shell® and the Binder Submalar®. Ideal for patients with both malar hypoplasia and midfacial soft-tissue deficiencies.

Cat. No.	Size	A	B	C	D
CSM-S	Small	4.38	3.00	0.40	0.30
CSM-M	Medium	4.90	3.26	0.50	0.40
CSM-L	Large	5.43	3.72	0.60	0.50
CSM-XL	X-Large	5.78	4.16	0.70	0.60
CSM-SIZ-NS	Sizes	Complete Set			

Dimensions in centimeters




Combined Malar



The Combined Malar Implants were designed to provide volume to the midfacial soft tissue with moderate anterior augmentation.

Catalog Number	Nominal Dimensions			
	A	B	C	D
S160-614S	4.5cm	3.0cm	0.4cm	0.35cm
S160-624S	5.0cm	3.3cm	0.45cm	0.4cm
S160-634S	5.6cm	3.7cm	0.5cm	0.45cm

- Sizer set S900-013
- Implant Supplied Sterile



SIMILAR PRODUCT


ALTERNATIVE



Conform™ Midfacial

Cat. No.	Size	A	B	C	D
CCSM-XS	X-Small	3.4	2.0	0.39	0.20
CCSM-S	Small	4.1	2.5	0.47	0.27
CCSM-M	Medium	4.7	2.7	0.53	0.30
CCSM-L	Large	5.1	3.0	0.56	0.32
CCSM-XL	X-Large	5.5	3.2	0.60	0.40
CCSM-SIZ-NS	Sizers	Complete Set			

Dimensions in centimeters




Anatomical Chin



Most similar to Implantech’s Extended Anatomical Chin with the XL providing more anterior projection.

Catalog Number	Nominal Dimensions				
	A	B	C	D	E
S250-114S	4.8cm	1.0cm	0.4cm	3.5cm	9.0cm
S250-124S	4.8cm	1.2cm	0.5cm	4.0cm	10.0cm
S250-134S	5.7cm	1.3cm	0.6cm	4.6cm	11.6cm
S250-144S	6.3cm	1.4cm	0.9cm	4.9cm	11.8cm

- Sizer set S900-001
- Implant Supplied Sterile



SIMILAR PRODUCT

ALTERNATIVE

ALTERNATIVE




Extended Anatomical Chin Implant

For additional 2mm anterior projection the GLE-2 wafer can be used

Cat. No.	Size	A	B	C	D	E	F
EAC-S	Small	5.5	1.0	0.45	9.0	3.7	0.2
EAC-M	Medium	5.6	1.2	0.50	9.5	4.2	0.4
EAC-L	Large	6.8	1.4	0.60	11.0	4.8	0.4
EAC-XL	X-Large	6.0	1.4	0.75	11.0	5.1	0.5
EAC-SIZ-NS	Sizers	Complete Set					

Dimensions in centimeters




Anatomical Chin



Most similar to Implantech’s Extended Anatomical Chin with the XL providing more anterior projection.

Catalog Number	Nominal Dimensions				
	A	B	C	D	E
S250-114S	4.8cm	1.0cm	0.4cm	3.5cm	9.0cm
S250-124S	4.8cm	1.2cm	0.5cm	4.0cm	10.0cm
S250-134S	5.7cm	1.3cm	0.6cm	4.6cm	11.6cm
S250-144S	6.3cm	1.4cm	0.9cm	4.9cm	11.8cm

- Sizer set S900-001
- Implant Supplied Sterile



SIMILAR PRODUCT

ALTERNATIVE


ALTERNATIVE



Conform™ Extended Anatomical Chin Implant

Cat. No.	Size	A	B	C	D	E	F
CEAC-S	Small	5.0	1.00	0.43	8.7	3.5	0.22
CEAC-M	Medium	5.6	1.20	0.50	10.5	4.2	0.25
CEAC-L	Large	6.2	1.35	0.60	11.5	4.5	0.31
CEAC-XL	X-Large	6.8	1.50	0.75	12.5	5.0	0.33
CEAC-SIZ-NS	Sizers	Complete Set					

Dimensions in centimeters






Most similar to Implantech’s Extended Anatomical Chin with the XL providing more anterior projection.

Catalog Number	Nominal Dimensions				
	A	B	C	D	E
S250-114S	4.8cm	1.0cm	0.4cm	3.5cm	9.0cm
S250-124S	4.8cm	1.2cm	0.5cm	4.0cm	10.0cm
S250-134S	5.7cm	1.3cm	0.6cm	4.6cm	11.6cm
S250-144S	6.3cm	1.4cm	0.9cm	4.9cm	11.8cm

- Sizer set S900-001
- Implant Supplied Sterile



SIMILAR PRODUCT

ALTERNATIVE


ALTERNATIVE




ePTFE-Coated Extended Anatomical Chin

Cat. No.	Size	A	B	C	D	E	F
ePTFE-EAC-S	Small	5.5	1.0	0.45	9.0	3.7	0.2
ePTFE-EAC-M	Medium	5.6	1.2	0.50	9.5	4.2	0.4
ePTFE-EAC-L	Large	6.8	1.4	0.60	11.0	4.8	0.4
ePTFE-EAC-XL	X-Large	6.0	1.4	0.75	11.0	5.1	0.5
EAC-SIZ-NS	Sizers	Complete Set (silicone)					

Dimensions in centimeters




Mittelman Geniomandibular Groove Chin/Jowl



Catalog Number	Nominal Dimensions				
	A	B	C	D	E
S250-264S	5.8cm	1.3cm	0.6cm	0.4cm	10.6cm
S250-274S	5.8cm	1.3cm	0.7cm	0.5cm	10.8cm
S250-284S	5.8cm	1.3cm	0.8cm	0.6cm	10.9cm

- Sizer set S900-023
- Implant Supplied Sterile




SIMILAR PRODUCT




Mandibular Pre Jowl
Chin™

Cat. No.	Size	A	B	C	D	E	F
MPJC-S	Small	6.0	1.5	0.65	10.5	4.0	0.33
MPJC-M	Medium	6.0	1.5	0.80	10.8	4.3	0.40
MPJC-L	Large	6.0	1.5	0.95	11.0	4.8	0.50
MPJC-XL	X-Large	6.0	1.5	1.10	11.0	5.0	0.56
MPJC-SIZ-NS	Sizers	Complete Set					

Dimensions in centimeters




Curvilinear Chin



Catalog Number	Nominal Dimensions		
	A	B	C
S250-314S	3.2cm	1.2cm	0.5cm
S250-324S	3.3cm	1.3cm	0.6cm
S250-334S	3.7cm	1.4cm	0.6cm
S250-344S	4.3cm	1.5cm	0.7cm
S250-354S	4.7cm	1.6cm	0.7cm
S250-364S	5.1cm	1.7cm	0.8cm

- Sizer set S900-007
- Implant Supplied Sterile



SIMILAR PRODUCT


ALTERNATIVE




Curvilinear Silicone Chin

Cat. No.	Size	A	B	C	D
CSC-XS	X-Small	3.1	1.1	0.45	3.0
CSC-S	Small	3.2	1.4	0.65	3.3
CSC-M	Medium	3.6	1.4	0.65	4.0
CSC-L	Large	4.4	1.4	0.65	4.8
CSC-XL	X-Large	4.0	1.5	0.70	4.5
CSC-SIZ-NS	Sizers	Complete Set			

Dimensions in centimeters




Curvilinear Chin



Catalog Number	Nominal Dimensions		
	A	B	C
S250-314S	3.2cm	1.2cm	0.5cm
S250-324S	3.3cm	1.3cm	0.6cm
S250-334S	3.7cm	1.4cm	0.6cm
S250-344S	4.3cm	1.5cm	0.7cm
S250-354S	4.7cm	1.6cm	0.7cm
S250-364S	5.1cm	1.7cm	0.8cm

- Sizer set S900-007
- Implant Supplied Sterile



SIMILAR PRODUCT


ALTERNATIVE



Anatomical Chin

Cat. No.	Size	A	B	C	D	E
AC-S	Small	4.7	1.2	0.60	6.0	2.1
AC-M	Medium	5.3	1.2	0.90	6.6	2.5
AC-L	Large	5.1	1.2	1.00	6.2	2.5
AC-XL	X-Large	5.3	1.2	1.20	6.6	2.8
AC-SIZ-NS	Sizers	Complete Set				

Dimensions in centimeters



Articulated Chin




250-414S, -424S

The Articulated Chin was designed to provide more anterior projection and also be more flexible due to the articulated back to reduce internal pressure on the mandible and provide a better fit.

Catalog Number	Nominal Dimensions				
	A	B	C	D	E
S250-414S	5.9cm	1.4cm	0.9cm	2.9cm	8.2cm
S250-424S	7.4cm	1.4cm	1.2cm	3.7cm	10.2cm

- Sizer set S900-004
- Implant Supplied Sterile



ALTERNATIVE

ALTERNATIVE

ALTERNATIVE


ALTERNATIVE



Flowers Mandibular Glove®

Cat. No.	Size	A	B	C	D	E	F
FMG-VXS	X-Small	6.00	1.75	0.50	9.0	3.65	0.2
FMG-VS	Small	6.02	1.75	0.60	9.0	3.75	0.3
FMG-VM	Medium	6.06	1.75	0.80	10.0	4.41	0.3
FMG-VL	Large	6.10	1.75	1.00	10.0	4.61	0.4
FMG-VXL	X-Large	6.30	1.75	1.40	10.8	5.40	0.5
FMG-VSIZ-NS	Sizers	Complete Set					

Dimensions in centimeters



If similar anterior projection is desired, consider using the Flowers Mandibular Glove® (FMG-L or FMG-XL) or the Mandibular Pre Jowl Chin™ (MPJC-XL). Made of silicone, these implants are already quite flexible.

If you desire more flexibility over more anterior projection, the Conform™ style EAC—with its gridlike backing—helps relax the implant's material memory, allowing it to better conform to the underlying tissues.

Articulated Chin




250-414S, -424S

The Articulated Chin was designed to provide more anterior projection and also be more flexible due to the articulated back to reduce internal pressure on the mandible and provide a better fit.

Catalog Number	Nominal Dimensions				
	A	B	C	D	E
S250-414S	5.9cm	1.4cm	0.9cm	2.9cm	8.2cm
S250-424S	7.4cm	1.4cm	1.2cm	3.7cm	10.2cm

- Sizer set S900-004
- Implant Supplied Sterile



ALTERNATIVE

ALTERNATIVE

ALTERNATIVE


ALTERNATIVE



Mandibular Pre Jowl
Chin™

Cat. No.	Size	A	B	C	D	E	F
MPJC-S	Small	6.0	1.5	0.65	10.5	4.0	0.33
MPJC-M	Medium	6.0	1.5	0.80	10.8	4.3	0.40
MPJC-L	Large	6.0	1.5	0.95	11.0	4.8	0.50
MPJC-XL	X-Large	6.0	1.5	1.10	11.0	5.0	0.56
MPJC-SIZ-NS	Sizers	Complete Set					

Dimensions in centimeters



If similar anterior projection is desired, consider using the Flowers Mandibular Glove® (FMG-L or FMG-XL) or the Mandibular Pre Jowl Chin™ (MPJC-XL). Made of silicone, these implants are already quite flexible.

If you desire more flexibility over more anterior projection, the Conform™ style EAC—with its gridlike backing—helps relax the implant’s material memory, allowing it to better conform to the underlying tissues.

Articulated Chin




250-414S, -424S

The Articulated Chin was designed to provide more anterior projection and also be more flexible due to the articulated back to reduce internal pressure on the mandible and provide a better fit.

Catalog Number	Nominal Dimensions				
	A	B	C	D	E
S250-414S	5.9cm	1.4cm	0.9cm	2.9cm	8.2cm
S250-424S	7.4cm	1.4cm	1.2cm	3.7cm	10.2cm

- Sizer set S900-004
- Implant Supplied Sterile



ALTERNATIVE

ALTERNATIVE

ALTERNATIVE


ALTERNATIVE



Extended Anatomical
Chin Implant

Cat. No.	Size	A	B	C	D	E	F
EAC-S	Small	5.5	1.0	0.45	9.0	3.7	0.2
EAC-M	Medium	5.6	1.2	0.50	9.5	4.2	0.4
EAC-L	Large	6.8	1.4	0.60	11.0	4.8	0.4
EAC-XL	X-Large	6.0	1.4	0.75	11.0	5.1	0.5
EAC-SIZ-NS	Sizers	Complete Set					

Dimensions in centimeters



If similar anterior projection is desired, consider using the Flowers Mandibular Glove® (FMG-L or FMG-XL) or the Mandibular Pre Jowl Chin™ (MPJC-XL). Made of silicone, these implants are already quite flexible.

If you desire more flexibility over more anterior projection, the Conform™ style EAC—with its gridlike backing—helps relax the implant’s material memory, allowing it to better conform to the underlying tissues.

Articulated Chin




250-414S, -424S

The Articulated Chin was designed to provide more anterior projection and also be more flexible due to the articulated back to reduce internal pressure on the mandible and provide a better fit.

Catalog Number	Nominal Dimensions				
	A	B	C	D	E
S250-414S	5.9cm	1.4cm	0.9cm	2.9cm	8.2cm
S250-424S	7.4cm	1.4cm	1.2cm	3.7cm	10.2cm

- Sizer set S900-004
- Implant Supplied Sterile



ALTERNATIVE

ALTERNATIVE

ALTERNATIVE


ALTERNATIVE



Conform™ Extended Anatomical Chin Implant

Cat. No.	Size	A	B	C	D	E	F
CEAC-S	Small	5.0	1.00	0.43	8.7	3.5	0.22
CEAC-M	Medium	5.6	1.20	0.50	10.5	4.2	0.25
CEAC-L	Large	6.2	1.35	0.60	11.5	4.5	0.31
CEAC-XL	X-Large	6.8	1.50	0.75	12.5	5.0	0.33
CEAC-SIZ-NS	Sizers	Complete Set					


Dimensions in centimeters



If similar anterior projection is desired, consider using the Flowers Mandibular Glove® (FMG-L or FMG-XL) or the Mandibular Pre Jowl Chin™ (MPJC-XL). Made of silicone, these implants are already quite flexible.

If you desire more flexibility over more anterior projection, the Conform™ style EAC—with its gridlike backing—helps relax the implant’s material memory, allowing it to better conform to the underlying tissues.


Symmetrical Chin



Another style of a chin implant to provide high central projection with extended "arms".

Catalog Number	Nominal Dimensions				
	A	B	C	D	E
S250-514S	5.7cm	1.3cm	0.6cm	2.3cm	8.0cm
S250-524S	6.3cm	1.5cm	0.8cm	3.9cm	9.9cm

- Sizer set S900-002
- Implant Supplied Sterile



ALTERNATIVE


ALTERNATIVE




Extended Anatomical
Chin Implant

Cat. No.	Size	A	B	C	D	E	F
EAC-S	Small	5.5	1.0	0.45	9.0	3.7	0.2
EAC-M	Medium	5.6	1.2	0.50	9.5	4.2	0.4
EAC-L	Large	6.8	1.4	0.60	11.0	4.8	0.4
EAC-XL	X-Large	6.0	1.4	0.75	11.0	5.1	0.5
EAC-SIZ-NS	Sizers	Complete Set					

Dimensions in centimeters




Symmetrical Chin



Another style of a chin implant to provide high central projection with extended "arms".

Catalog Number	Nominal Dimensions				
	A	B	C	D	E
S250-514S	5.7cm	1.3cm	0.6cm	2.3cm	8.0cm
S250-524S	6.3cm	1.5cm	0.8cm	3.9cm	9.9cm

- Sizer set S900-002
- Implant Supplied Sterile



ALTERNATIVE


ALTERNATIVE



Mandibular Pre Jowl
Chin™

Cat. No.	Size	A	B	C	D	E	F
MPJC-S	Small	6.0	1.5	0.65	10.5	4.0	0.33
MPJC-M	Medium	6.0	1.5	0.80	10.8	4.3	0.40
MPJC-L	Large	6.0	1.5	0.95	11.0	4.8	0.50
MPJC-XL	X-Large	6.0	1.5	1.10	11.0	5.0	0.56
MPJC-SIZ-NS	Sizers	Complete Set					

Dimensions in centimeters




Cleft Chin



Catalog Number	Nominal Dimensions					
	A	B	C	D	E	F
S250-614S	4.7cm	1.3cm	0.7cm	3.2cm	8.2cm	0.5cm
S250-624S	5.6cm	1.3cm	0.9cm	4.1cm	9.9cm	0.7cm
S250-634S	7.0cm	1.3cm	1.0cm	4.6cm	11.4cm	0.8cm

- Sizer set S900-005
- Implant Supplied Sterile



ALTERNATIVE


ALTERNATIVE



Extended Anatomical
Chin Implant

Cat. No.	Size	A	B	C	D	E	F
EAC-S	Small	5.5	1.0	0.45	9.0	3.7	0.2
EAC-M	Medium	5.6	1.2	0.50	9.5	4.2	0.4
EAC-L	Large	6.8	1.4	0.60	11.0	4.8	0.4
EAC-XL	X-Large	6.0	1.4	0.75	11.0	5.1	0.5
EAC-SIZ-NS	Sizers	Complete Set					

Dimensions in centimeters




Important note: Chin clefts are created by a divit in the soft tissue and not bone structure. The EAC or FMG can be notched and then soft tissue manipulation done to create the cleft.



Catalog Number	Nominal Dimensions					
	A	B	C	D	E	F
S250-614S	4.7cm	1.3cm	0.7cm	3.2cm	8.2cm	0.5cm
S250-624S	5.6cm	1.3cm	0.9cm	4.1cm	9.9cm	0.7cm
S250-634S	7.0cm	1.3cm	1.0cm	4.6cm	11.4cm	0.8cm


- Sizer set S900-005
- Implant Supplied Sterile



Flowers Mandibular Glove®

Cat. No.	Size	A	B	C	D	E	F
FMG-VXS	X-Small	6.00	1.75	0.50	9.0	3.65	0.2
FMG-VS	Small	6.02	1.75	0.60	9.0	3.75	0.3
FMG-VM	Medium	6.06	1.75	0.80	10.0	4.41	0.3
FMG-VL	Large	6.10	1.75	1.00	10.0	4.61	0.4
FMG-VXL	X-Large	6.30	1.75	1.40	10.8	5.40	0.5
FMG-VSIZ-NS	Sizers	Complete Set					

Dimensions in centimeters




Important note: Chin clefts are created by a divit in the soft tissue and not bone structure. The EAC or FMG can be notched and then soft tissue manipulation done to create the cleft.

Mittelman Geniomandibular Groove Jowl



Catalog Number	Nominal Dimensions				
	A	B	C	D	E
S250-764S	6.0cm	1.4cm	0.1cm	0.4cm	10.5cm
S250-774S	6.0cm	1.4cm	0.1cm	0.5cm	10.5cm
S250-784S	6.0cm	1.4cm	0.1cm	0.6cm	10.5cm

- Sizer set S900-024
- Implant Supplied Sterile



SIMILAR PRODUCT




Mandibular Pre Jowl®

Provides anterior chin augmentation with a subtle increase in lateral fullness to augment the pre-jowl sulcus.

Cat. No.	Size	A	B	C	D	E	F
MPJ-S	Small	6.0	1.0	0.10	10.9	4.2	0.3
MPJ-M	Medium	6.0	1.0	0.10	11.0	4.3	0.4
MPJ-L	Large	6.0	1.0	0.10	11.2	4.3	0.5
MPJ-XL	X-Large	6.0	1.0	0.10	12.5	4.6	0.6
MPJ-SIZ-NS	Sizers	Complete Set					

Dimensions in centimeters




Bilateral Groove Chin



Catalog Number	Nominal Dimensions					
	A	B	C	D	E	F
S250-814S	5.8cm	1.6cm	0.2cm	3.5cm	9.6cm	0.5cm
S250-824S	5.8cm	1.6cm	0.4cm	3.7cm	9.7cm	0.5cm
S250-834S	5.9cm	1.6cm	0.5cm	4.3cm	10.6cm	0.8cm
S250-844S	5.9cm	1.6cm	0.7cm	4.4cm	10.8cm	1.0cm
S250-854S	6.1cm	1.7cm	1.2cm	5.3cm	12.0cm	1.5cm

- Sizer set S900-003
- Implant Supplied Sterile




SIMILAR PRODUCT




Flowers Mandibular Glove®

Cat. No.	Size	A	B	C	D	E	F
FMG-VXS	X-Small	6.00	1.75	0.50	9.0	3.65	0.2
FMG-VS	Small	6.02	1.75	0.60	9.0	3.75	0.3
FMG-VM	Medium	6.06	1.75	0.80	10.0	4.41	0.3
FMG-VL	Large	6.10	1.75	1.00	10.0	4.61	0.4
FMG-VXL	X-Large	6.30	1.75	1.40	10.8	5.40	0.5
FMG-VSIZ-NS	Sizes	Complete Set					

Dimensions in centimeters




Concave Back Chin



Catalog Number	Nominal Dimensions		
	A	B	C
S250-914S	3.4cm	1.3cm	0.55cm
S250-924S	3.6cm	1.4cm	0.6cm
S250-934S	4.1cm	1.5cm	0.7cm
S250-944S	4.7cm	1.7cm	1.0cm

- Sizer set S900-006
- Implant Supplied Sterile



SIMILAR PRODUCT




Curvilinear Silicone Chin

The biggest difference is that Implantech's CSC does not have varying height to the implant.

Cat. No.	Size	A	B	C	D
CSC-XS	X-Small	3.1	1.1	0.45	3.0
CSC-S	Small	3.2	1.4	0.65	3.3
CSC-M	Medium	3.6	1.4	0.65	4.0
CSC-L	Large	4.4	1.4	0.65	4.8
CSC-XL	X-Large	4.0	1.5	0.70	4.5
CSC-SIZ-NS	Sizers	Complete Set			

Dimensions in centimeters




Spectrum Dorsal Columella Styles 1-4



Catalog Number	Nominal Dimensions					
	A	B	C	D	E	F
S400-114S	5.8cm	3.3cm	0.6cm	0.7cm	0.4cm	0.5cm
S400-214S	5.7cm	2.5cm	0.8cm	0.6cm	0.5cm	0.3cm
S400-314S	5.2cm	3.0cm	0.9cm	0.8cm	0.5cm	0.3cm
S400-414S	5.0cm	2.0cm	0.9cm	0.6cm	0.4cm	0.4cm

- Sizer set(s) S900-016, S900-020, S900-021, S900-022
- Implant Supplied Sterile



ALTERNATIVE

ALTERNATIVE




Voloshin Dorsal Columella

Voloshin Dorsal Columella (VSDC) or Dorsal Columella (DCI). Suggest DCI for longer lengths. This is a dorsal columella template meant to be carved into the specific shape needed. Available in two different durometers: Firm and Soft

Cat. No.	A	B	C	D	E	F	G	H
VSDC-001	5.0	3.0	1.1	0.7	0.4	0.22	0.40	0.4
VSDC-002	5.4	3.0	1.0	0.7	0.4	0.20	0.40	0.4
VSDC-003	4.9	3.2	1.0	0.7	0.4	0.43	0.70	0.6
VSDC-004	5.2	3.2	1.2	0.7	0.4	0.50	0.75	0.6
VSDC-010	5.2	3.3	1.0	0.7	0.4	0.35	0.55	0.6
VSDC-SIZ-NS	Sizers	Complete Set						

Dimensions in centimeters




Spectrum Dorsal Columella Styles 1-4



Catalog Number	Nominal Dimensions					
	A	B	C	D	E	F
S400-114S	5.8cm	3.3cm	0.6cm	0.7cm	0.4cm	0.5cm
S400-214S	5.7cm	2.5cm	0.8cm	0.6cm	0.5cm	0.3cm
S400-314S	5.2cm	3.0cm	0.9cm	0.8cm	0.5cm	0.3cm
S400-414S	5.0cm	2.0cm	0.9cm	0.6cm	0.4cm	0.4cm

- Sizer set(s) S900-016, S900-020, S900-021, S900-022
- Implant Supplied Sterile



ALTERNATIVE




Dorsal Columella

Voloshin Dorsal Columella (VSDC) or Dorsal Columella (DCI). Suggest DCI for longer lengths. This is a dorsal columella template meant to be carved into the specific shape needed. Available in two different durometers: Firm and Soft

Cat. No.	Type	A	B	C	D	E	F	G
DCI-03F	Firm	5.8	3.0	1.1	0.3	0.9	0.8	0.6
DCI-03S	Soft	5.8	3.0	1.1	0.3	0.9	0.8	0.6

Dimensions in centimeters




Spectrum Dorsal Columella Style 5



Catalog Number	Nominal Dimensions					
	A	B	C	D	E	F
S400-514S	4.8cm	2.9cm	0.9cm	0.6cm	0.5cm	0.4cm
S400-524S	4.5cm	2.9cm	1.0cm	0.6cm	0.4cm	0.4cm
S400-534S	4.8cm	2.9cm	1.1cm	0.6cm	0.4cm	0.4cm
S400-544S	4.8cm	2.9cm	1.2cm	0.6cm	0.5cm	0.5cm

- Sizer set S900-017
- Implant Supplied Sterile



ALTERNATIVE


ALTERNATIVE



Voloshin Dorsal Columella

Cat. No.	A	B	C	D	E	F	G	H
VSDC-001	5.0	3.0	1.1	0.7	0.4	0.22	0.40	0.4
VSDC-002	5.4	3.0	1.0	0.7	0.4	0.20	0.40	0.4
VSDC-003	4.9	3.2	1.0	0.7	0.4	0.43	0.70	0.6
VSDC-004	5.2	3.2	1.2	0.7	0.4	0.50	0.75	0.6
VSDC-010	5.2	3.3	1.0	0.7	0.4	0.35	0.55	0.6
VSDC-SIZ-NS	Sizers	Complete Set						


Dimensions in centimeters





Catalog Number	Nominal Dimensions					
	A	B	C	D	E	F
S400-514S	4.8cm	2.9cm	0.9cm	0.6cm	0.5cm	0.4cm
S400-524S	4.5cm	2.9cm	1.0cm	0.6cm	0.4cm	0.4cm
S400-534S	4.8cm	2.9cm	1.1cm	0.6cm	0.4cm	0.4cm
S400-544S	4.8cm	2.9cm	1.2cm	0.6cm	0.5cm	0.5cm


- Sizer set S900-017
- Implant Supplied Sterile



Anatomical Nasal Implant

Cat. No.	A	B	C	D	E	F	G	H
ANI-010	3.7	2.3	0.80	0.66	0.80	0.20	0.58	0.42
ANI-020	4.0	2.5	0.90	0.72	0.90	0.21	0.65	0.45
ANI-030	4.4	2.7	1.00	0.79	0.95	0.23	0.70	0.50
ANI-031	4.4	2.7	1.00	0.79	0.95	0.42	0.70	0.50
ANI-040	4.7	2.9	1.05	0.84	1.00	0.23	0.75	0.55
ANI-041	4.7	2.9	1.05	0.84	1.00	0.45	0.75	0.55
ANI-SIZ-NS	Sizers	Complete Set						


Dimensions in centimeters





Catalog Number	Nominal Dimensions					
	A	B	C	D	E	F
S400-614S	4.7cm	2.9cm	1.1cm	0.6cm	0.3cm	0.4cm
S400-624S	4.7cm	3.0cm	1.1cm	0.5cm	0.5cm	0.4cm
S400-634S	4.8cm	3.0cm	1.1cm	0.7cm	0.5cm	0.4cm
S400-644S	4.8cm	3.0cm	1.2cm	0.6cm	0.6cm	0.4cm
S400-654S	4.6cm	3.0cm	1.0cm	0.7cm	0.4cm	0.4cm
S400-664S	5.0cm	3.0cm	1.1cm	0.7cm	0.3cm	0.4cm


- Sizer set S900-015
- Implant Supplied Sterile



Voloshin Dorsal Columella

Cat. No.	A	B	C	D	E	F	G	H
VSDG-001	5.0	3.0	1.1	0.7	0.4	0.22	0.40	0.4
VSDG-002	5.4	3.0	1.0	0.7	0.4	0.20	0.40	0.4
VSDG-003	4.9	3.2	1.0	0.7	0.4	0.43	0.70	0.6
VSDG-004	5.2	3.2	1.2	0.7	0.4	0.50	0.75	0.6
VSDG-010	5.2	3.3	1.0	0.7	0.4	0.35	0.55	0.6
VSDG-SIZ-NS	Sizes	Complete Set						


Dimensions in centimeters





Catalog Number	Nominal Dimensions					
	A	B	C	D	E	F
S400-614S	4.7cm	2.9cm	1.1cm	0.6cm	0.3cm	0.4cm
S400-624S	4.7cm	3.0cm	1.1cm	0.5cm	0.5cm	0.4cm
S400-634S	4.8cm	3.0cm	1.1cm	0.7cm	0.5cm	0.4cm
S400-644S	4.8cm	3.0cm	1.2cm	0.6cm	0.6cm	0.4cm
S400-654S	4.6cm	3.0cm	1.0cm	0.7cm	0.4cm	0.4cm
S400-664S	5.0cm	3.0cm	1.1cm	0.7cm	0.3cm	0.4cm


- Sizer set S900-015
- Implant Supplied Sterile




ePTFE-Coated Voloshin
Dorsal Columella

Cat. No.	A	B	C	D	E	F	G	H
ePTFE-VSDG-001	5.0	3.0	1.1	0.7	0.4	0.22	0.40	0.4
ePTFE-VSDG-002	5.4	3.0	1.0	0.7	0.4	0.20	0.40	0.4
ePTFE-VSDG-003	4.9	3.2	1.0	0.7	0.4	0.43	0.70	0.6
ePTFE-VSDG-004	5.2	3.2	1.2	0.7	0.4	0.50	0.75	0.6
ePTFE-VSDG-010	5.2	3.3	1.0	0.7	0.4	0.35	0.55	0.6
VSDG-SIZ-NS	Sizers	Complete Set (silicone)						

Dimensions in centimeters




Nasal Dorsum Style 7



Catalog Number	Nominal Dimensions			
	A	B	C	D
S400-714S	2.8cm	0.7cm	0.3cm	0.4cm
S400-724S	3.2cm	0.8cm	0.3cm	0.4cm
S400-734S	3.8cm	0.9cm	0.4cm	0.5cm
S400-744S	4.7cm	0.9cm	0.4cm	0.5cm
S400-754S	5.4cm	0.9cm	0.4cm	0.5cm

- Sizer set S900-018
- Implant Supplied Sterile



SIMILAR PRODUCT

ALTERNATIVE




Flowers Dorsal Nasal


Similar to the medium length of the Nasal Dorsum Style 7. If a shorter length is needed, simply trim off the excess.

Cat. No.	Size	A	B	C
FDI-1.5	Small	4.3	1.00	0.15
FDI-2.5	Sm-Med	4.4	1.05	0.25
FDI-3.5	Medium	4.4	1.10	0.35
FDI-4.5	Large	4.4	1.20	0.45
FDI-SIZ-NS	Sizers	Complete Set		

Dimensions in centimeters




Nasal Dorsum Style 7



Catalog Number	Nominal Dimensions			
	A	B	C	D
S400-714S	2.8cm	0.7cm	0.3cm	0.4cm
S400-724S	3.2cm	0.8cm	0.3cm	0.4cm
S400-734S	3.8cm	0.9cm	0.4cm	0.5cm
S400-744S	4.7cm	0.9cm	0.4cm	0.5cm
S400-754S	5.4cm	0.9cm	0.4cm	0.5cm

- Sizer set S900-018
- Implant Supplied Sterile



SIMILAR PRODUCT

ALTERNATIVE



ePTFE-Coated Flowers
Dorsal Nasal

Similar to the medium length of the Nasal Dorsum Style 7. If a shorter length is needed, simply trim off the excess.

Cat. No.	Size	A	B	C
ePTFE-FDI-1.5	Small	4.3	1.00	0.15
ePTFE-FDI-2.5	Sm-Med	4.4	1.05	0.25
ePTFE-FDI-3.5	Medium	4.4	1.10	0.35
ePTFE-FDI-4.5	Large	4.4	1.10	0.45
FDI-SIZ-NS	Sizers	Complete Set (silicone)		

Dimensions in centimeters

